



# EndoRing® II GelWell® Cups

## Single Use GelWell Cups

The single use GelWell cups are used for storing lubricating gels and liquids while using the EndoRing II. The single use cups are secured by an adapter ring optionally attached to the EndoRing foam cup. Single use convenience allows the operator to simply dispose of the cup when the procedure is complete.

### TECHNICAL SUMMARY

**EndoRing II:** Assemblies are molded from an autoclavable blue thermoplastic resin. Assembly consists of six pieces: 1 EndoRing Cup, 3 Finger Rings (S, M, L), 1 reusable GelWell Cup, and 1 GelWell Adapter Ring. Note: the EndoRing Cup has a molded plastic 30-mm scale for right or left handed users.

**EndoRing Metal Ruler:** Autoclavable anodized blue aluminum detachable ruler with contrasting white laser-etched scales in 1mm and 1/2mm increments. Accessory scale provides marking up to 37mm for extra long canals. Ruler can be used by right and left handed users. Available with or without rubber stop locks.

**GelWell System:** Reusable Cup: 2ml cup is molded from an autoclavable blue thermoplastic resin and attaches to the EndoRing II organizer. Single Use Cups: 0.65ml cups molded from a clear thermoplastic resin. The single use cup is placed into the holder and attaches to the EndoRing II organizer.


**e-Foam:** Single use patented open cell yellow and blue foam. Can be steam autoclave processed prior to use. During steam autoclave processing, endodontic instruments can be inserted into the foam.

### EndoRing II Endodontic Instrument Family of Products

The EndoRing II organizer enables the clinician to directly place, store, measure and clean endodontic hand and rotary instruments within the operating field. It reduces the risk of cross-contamination and minimizes the passing of sharps between dentists and staff. Using the EndoRing II organizer helps make endodontic procedures safer and faster.



Discover our complete line of quality Jordco® products at [www.jordco.com](http://www.jordco.com) or contact your dental supply dealer.

 Jordco, Incorporated  
8531 NE Cornell Rd, STE 800, Hillsboro, Oregon 97124 USA  
(800) 752-2812 toll free USA [www.jordco.com](http://www.jordco.com)  
Visit our website for patent and trademark information.



### PRODUCT SPECS

**Name:** EndoRing II GelWell Single Use Cups  
**REF:** ERG-s  
**Weight:** 2 oz

### PACKAGE CONTENTS

- 48 Single Use GelWell Cups

### FEATURES AND BENEFITS

#### Hand-held Instrument

Allows direct operator control of hand and rotary instruments. Eliminates time consuming passing of dangerous files from assistant or bracket tray. Minimizes eyestrain and loss of attention from searching for files, etc. Saves several minutes each procedure.

#### Two Color EndoFoam Cleans and Organizes

Jordco blue and yellow e-Foam inserts provide an additional option for color coding instrument set-ups. While in use, e-Foam efficiently cleans bioburden from most all endodontic files.

#### Detachable Laser Etched Metal Ruler

Autoclavable metal ruler with a precise laser etched right or left-handed 30mm scale. Includes 1/2 mm calibration markings on the lower half of the scale. The metal ruler also includes a secondary scale to measure files up to 37-mm in length. Includes a calibration track for easy positioning of instruments. Ruler can be easily detached for cleaning and processing.

#### Easy Calibration of Curved Instruments

Curved files can be easily positioned into the calibrating track within the standard or metal ruler and temporarily straightened with slight pressure. After the stop is adjusted, the file is removed and allowed to reassume its curved shape.

#### Single Use GelWell cups

GelWell cups offer the convenience of holding all lubricating gels including Jordco's EndoGel. The unique conical design allows the operator to easily load any endodontic file with EndoGel. With single use convenience, it makes clean up easy and helps to reduce cross contamination.